Rollautoklaven Rotating autoclaves

1231

Allgemeines

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General

Rotating autoclaves give a good mixing of the reaction components i. e. liquid-catalyst-gas. They are used in all cases where autoclaves with stirrer units type 1220, 1221, 1222, are not possible because chemical processes or design do not allow it.


Design rotating autoclaves

Autoclaves are supplied with a screw closure or lid. Standard fittings are gas inlet/outlet valve, pressure gauge, thermometer tube or thermometer bore and bursting disc. The autoclaves are electrically heated with heating elements built in a detachable housing.

The autoclave flange runs on rollers, and the body is slid into a coupling at the base. The driving shaft is well supported and driven by a belt from a geared electric motor to give approximately 40-70 rpm. It is possible to incline the autoclave by hand or hydraulically.

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Rotating autoclave example:
0.5 ltrs., 325 bar, 300 °C
cladged with tantalum

Material

Standard autoclaves are produced from stainless steel (Cr-Ni-Mo-steel), rust and acid resistant. Special constructions are possible in other materials in part by cladding e. g. ferritic steel, copper, nickel, monel, inconel, hastelloy, tantalum, titanium, zirconium, rhenium, silver, gold, palladium, platinum etc.

Pressure Range

Autoclaves can be used in positive pressure and also in vacuum ranges. Standard is: 50, 100, 200, 325, 400 and 500 bar. Variations are possible. The pressure stages are in accordance with DIN 2401, based on the guiding principals of the chemical industry **10/20/50/100/160/200/250/325/400/500/600/700/1000/1600/2500** bar. The marked values are the most common.

Temperature Range

The desired material of construction defines the temperature limitation. It generally holds that pressure and temperature must be balanced against one another. Temperatures below zero are possible mainly as special construction. Standard temperature is 350 °C and 500 °C.


Inspection

The inspection rules of the autoclaves are fixed in German Pressure Vessel-code, in which specified tests for material measurement and pressure tests by the German Technical Supervisory Association are considered.

Autoclaves for export can be certified for other acceptance organisations e. g. Dienst van het Stoomwezen, Büro Veritas, British Lloyds Register of Shipping, Shwedisch SA and others.

Furthermore autoclaves delivered with electric heating will be pressure and temperature checked (max. operating conditions) for tightness and faults.

Technical alterations reserved.

Copyright by Ernst Haage, Mülheim (Ruhr) 11/96

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Estanit GmbH · Zinkhüttenstraße 17 · D-45473 Mülheim an der Ruhr