Schüttelautoklaven Shaking autoclaves

1230

Allgemeines

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General

Shaking autoclaves give a good mixing of the reaction components i. e. liquid-catalyst-gas. They are used in all cases where autoclaves with stirrer units type 1220, 1221, 1222, are not possible because chemical processes or design do not allow it.

Design


Autoclaves are supplied with a screw closure or lid. Standard fittings are gas inlet/outlet valve, pressure gauge, thermometer tube or thermometer bore and bursting disc. A covibrative capillary for gas inlet is possible. A reflux condenser (on top of the lid) for gas sampling during reaction can be provided.

Shaking autoclaves are normally electrically heated. They are inserted in the heater housing and then fastened. Jacket heating with steam or oil is difficult because of the shaking motion. This is therefore not recommended. Electrical heating in explosion-proof design is possible.

For driving, a geared electric motor is used to give approximately 40-50 rpm. The rotative movement is transmitted to the autoclave by an eccentric. The autoclave is supported at the centre of gravity in a fork. This arrangement allows a biaxial movement which guarantees good mixing.

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Material

Standard autoclaves are produced from stainless steel (Cr-Ni-Mo-steel), rust and acid resistant. Special constructions are possible in other materials in part by cladding e. g. ferritic steel, copper, nickel, monel, inconel, hastelloy, tantalum, titanium, zirconium, rhenium, silver, gold, palladium, platinum etc.

Pressure Range

Autoclaves can be used in positive pressure and also in vacuum ranges. Standard is: 50, 100, 200, 325, 400 and 500 bar. Variations are possible. The pressure stages are in accordance with DIN 2401, based on the guiding principals of the chemical industry **10/20/50/100/160/200/250/325/400/500/600/700/1000/1600/2500** bar. The marked values are the most common.

Temperature Range

The desired material of construction defines the temperature limitation. It generally holds that pressure and temperature must be balanced against one another. Temperatures below zero are possible mainly as special construction. Standard temperature is 350 °C and 500 °C.


Inspection

The inspection rules of the autoclaves are fixed in German Pressure Vessel-code, in which specified tests for material measurement and pressure tests by the German Technical Supervisory Association are considered.

Autoclaves for export can be certified for other acceptance organisations e. g. Dienst van het Stoomwezen, Büro Veritas, British Lloyds Register of Shipping, Swedish SA and others.

Furthermore autoclaves delivered with electric heating will be pressure and temperature checked (max. operating conditions) for tightness and faults.

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Shaking autoclave example:
2 ltrs., 120 bar, 300 °C
cladded with silver

Copyright by Ernst Haage, Mülheim (Ruhr) 11/96

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Estanit GmbH · Zinkhüttenstraße 17 · D-45473 Mülheim an der Ruhr